

CAPI • DRAW

CAPILLARY BLOOD COLLECTION SYSTEM

- **NO SCOOPING** reduces hemolysis
- **ANTICOAGULANT COATED CAPILLARY AND MICRO-TUBE**—provides immediate mixing to prevent clotting
- **PRE-ATTACHED CAP WITH RE-SEALABLE SAMPLE PORT**—for ease of use and safety
- **ACCURATE SAMPLE VOLUME**—provides the proper sample / additive ratio
- **UNBREAKABLE, FLAT-BOTTOM MICRO-TUBE**—also fits most commonly available centrifuge adapters.



THE MOST COMPLETE RANGE OF ADDITIVES

(Custom configurations available)



Single-Use Only • In-Vitro Diagnostic
Non-Sterile • Latex-Free • 100% Plastic • Ready to Use

50 tubes / bag 10 bags / case (500 tubes)

CHEMISTRY

- 505 SERUM SEPARATOR GEL 200uL (red)
- 515 SERUM BILIRUBIN GEL 200uL (red with amber vial)
- 525 LITHIUM HEPARIN 200uL (green)
- 535 LITHIUM HEPARIN GEL 200uL (green)
- 545 SODIUM FLUORIDE 200uL (yellow)

HEMATOLOGY

- 510 EDTA 200uL (lavender)
- 511 EDTA 125uL (lavender)

COAGULATION

- 550 SODIUM CITRATE 3.2% 250uL (blue)



INNOVATIVE
Med Tech

866-560-1820 fax: 800-768-2825

www.innovativemedtech.com



Cross Reference

CAPI-DRAW CAPILLARY COLLECTOR			
Vendor	Vendor #	Description	CAP COLOR
IMT	511	Tube Capi-Draw EDTA 125ul 50/ Bx	LAV
RAM	07 6011	TUBE CAP COLL W/EDTA 125UL	LAV
STATSPIN	SS1E	STATSAMPLER 100UL EDTA	LAV
IMT	510	Tube Capi-Draw EDTA 200ul 50/ Bx	LAV
RAM	07 7051	TUBE CAP COLL EDTA 200UL	LAV
STATSPIN	SS2E	STATSAMPLER 200UL EDTA	LAV
SAR	20.1288.100	MICROTAINER MICROVETTE .2ML	LAV
IMT	505	Tube Capi-Draw Serum Gel 50/ Bx	RED
RAM	07 7120	TUBE CAP COLL SERUM GEL200UL	RED
STATSPIN	SS2U	STATSAMPLER SERUM GEL 200UL	RED
SAR	20.1290.100	MICROVETTE CB200 SERUM 100BX	RED
IMT	525	Tube Capi Draw Lith Hep 50/ Bx	GRN
RAM	07 7250	TUBE CAP W/HEP 200UL 50/BX	GRN
STATSPIN	SS2H	STATSAMPLER 200UL LI HEP	GRN
SAR	20.1292.100	MICROVETTE CB200 LI HEP 100B	GRN

GRAVITY-FLOW COLLECTORS			
Vendor	Vendor #	Description	CAP COLOR
IMT	510	Tube Capi-Draw EDTA 200ul 50/ Bx	LAV
BDV	365973	TUBE MICROTAINER W/EDTA	LAV
BDV	365974	TUBE MICROTAINER LAV W/EDTA	LAV
TER	T-MQK	TUBE CAPIJECT BLOOD LAVENDR	LAV
IMT	505	Tube Capi-Draw Serum Gel 50/ Bx	RED
BDV	365956	TUBE MICROTAINER SST	RED
BDV	365963	TUBE MICROTAINER RED 50/BX	RED
BDV	365967	TUBE MICROTAINER SST CLEAR	RED
TER	T-MG	TUBE CAPIJECT BLOOD CHRY RED	RED
TER	T-MGA	TUBE CAPIJECT BLOOD GOLD	RED
IMT	3095	Tube Capi Draw Lith Hep 50/ Bx	GRN
BDV	365971	TUBE MICROTAINER LITH HEP	GRN
BDV	365958	TUBE MICROTAINER W/LITH HEP	GRN
TER	PER-T-MLH	TUBE CAPIJECT GRN 1000/C	GRN