

BLOOD-GAS SAMPLING SYRINGE

VETERINARY AND LAB USE ONLY



Blood Gas Sampling Syringes serve as arterial extraction devices for pH/blood gas, electrolyte, oximetry and metabolite analysis. The devices are coated with approx. 50 I.U. balanced heparin. Sterile, single use. Compatible with all common blood-gas analyzer types.

<u>REF</u>	<u>Product description</u>	<u>Box</u>	<u>Case</u>
3999599	2 ml, Luer Slip, sterile, ~50 I. U. bal Heparin	50 ea	1000 ea
999017	1 ml, Luer Slip, sterile, ~50 I. U. bal Heparin	50 ea	1000 ea

INNOVATIVE
Med Tech 

866-560-1820
www.InnovativeMedTech.com